

This disinfectant is a dual chain quaternary ammonium/high level alcohol solution. Use in Health Care & Alternate Care Settings such as hospitals, intensive care units, emergency medical settings, central supply, laboratories, clinics, nursing homes, and dental suites where control of the cross-contamination hazards between treated surfaces is required. It is effective against the following pathogens on hard, non-porous contact surfaces when used as directed:  
 Hepatitis B Virus\*\* (HBV) (1 min.), Hepatitis C Virus\*\* (HCV) (1 min.), Human Coronavirus\*\* (1 min.), HIV-1\*\* (AIDS virus) (1 min.), Staphylococcus aureus (Staph) (1 min.), Pseudomonas aeruginosa (1 min.), Salmonella choleraesuis (1 min.), Escherichia coli (E. coli) (1 min.), Adenovirus type 2\*\* (1 min.), Methicillin Resistant Staphylococcus aureus (MRSA) (1 min.), Vancomycin Resistant Enterococcus faecalis (VRE) (1 min.), Trichophyton mentagrophytes (Athlete's Foot fungus) (1 min.), Herpes Simplex Virus type 2\*\* (genital herpes virus) (1 min.), Mycobacterium bovis (Tuberculosis bacteria) (1 min. at 20°C)

\* A precleaning step is required for HIV-1, HBV and HCV.  
 \*\* A precleaning step is required.

**DIRECTIONS FOR USE**

It is a violation of Federal law to use this product in a manner inconsistent with its labeling.

**GENERAL PRECAUTIONS AND RESTRICTIONS:** This product is not to be used as a terminal sterilant/high level disinfectant on any surface or instrument that (1) is introduced directly into the human body, either into or in contact with the bloodstream or normally sterile areas of the body, or (2) contact intact mucous membranes but which does not ordinarily penetrate the blood barrier or otherwise enter normally sterile areas of the body. This product may be used to preclean or decontaminate critical or semi-critical medical devices prior to sterilization or high level disinfection.

**TO CLEAN, DECONTAMINATE, AND DISINFECT:** Thoroughly wet the surface with this product and allow to remain wet for 1 minute. When used as a fungicidal disinfectant, preclean the surfaces to be disinfected. This product can be used for this purpose. Then, thoroughly wet the surface with this product and allow to remain wet for 1 minute.

**For Disinfection of:** inanimate, hard, non-porous surfaces such as chairs, stools, counter tops, drawer pulls, door knobs and handles, carts, baskets, tables, cabinets, lighting fixtures, bassinets, bed rails, incubators, physical therapy equipment, shower stalls, toilet exteriors and toilet seats, telephones and other related hard surfaces in hospitals, operating and surgical rooms, intensive care and critical care units, emergency units, dental and veterinary offices, surgery, recovery, anesthesia, X-Ray, Cat. Lab., blood collection and donation centers, isolation units, orthopedics, newborn nursery, respiratory therapy, radiology, pediatrics, physical therapy, rehabilitation, ophthalmology, optometry, hospices, central supply, emergency medical vehicles, police and fire vehicles, correctional facilities, funeral homes, morgues, spas, health clubs, salons, schools and school medical facilities, school buses, day care and non-food contact surfaces.

Use on stainless steel, glass, formica and other similar surfaces. Use care when using on varnished surfaces or plexi-glass.



Kills HIV-1, HBV, and HCV on PRE-CLEANED ENVIRONMENTAL SURFACES/ OBJECTS PREVIOUSLY SOILED WITH BLOOD/BODY FLUIDS in health care settings (Hospitals, Nursing Homes) or other settings in which there is an expected likelihood of soiling of inanimate surfaces/objects with blood or body fluids, and in which the surfaces/objects likely to be soiled with blood or body fluids can be associated with the potential for transmission of Human Immunodeficiency Virus Type 1 (HIV-1) (associated with AIDS), Hepatitis B Virus (HBV), and Hepatitis C Virus (HCV).

**SPECIAL INSTRUCTIONS FOR CLEANING AND DECONTAMINATION AGAINST HIV-1 (AIDS VIRUS), HBV, and HCV FOR SURFACES OR OBJECTS SOILED WITH BLOOD OR BODY FLUIDS.**

**PERSONAL PROTECTION:** Wear disposable latex gloves, protective gown, face mask and eye covering, as appropriate, when handling HIV-1, HBV, and HCV infected blood or body fluids. **CLEANING PROCEDURE:** All blood and other body fluids must be thoroughly cleaned from surfaces and objects before application of this disinfectant. This product can be used for this purpose. **CONTACT TIME:** Thoroughly wet surface with this

product. Allow to air dry. The efficacy of a 1 minute contact time has been shown to be adequate against HIV-1 (AIDS virus), HBV, and HCV.

**DISPOSAL OF INFECTIOUS MATERIALS:** Blood and other body fluids should be autoclaved and disposed of according to Federal, State, and local regulations for infectious waste disposal.

**Precautionary Statements:** Hazards to Humans and Domestic Animals: **WARNING.** Causes substantial but temporary eye injury. Do not get in eyes or on clothing. Wear protective eyewear. Wash thoroughly with soap and water after handling. Remove contaminated clothing and wash clothing before reuse.

**Physical or Chemical Hazards:** Flammable. Keep away from heat or open flame.

First Aid	
<b>If in eyes:</b>	<ul style="list-style-type: none"> <li>• Hold eye open and rinse slowly and gently with water for 15-20 minutes.</li> <li>• Remove contact lenses, if present, after the first 5 minutes, then continue rinsing eye.</li> <li>• Call a poison control center or doctor for treatment advice.</li> </ul>
Have the product container or label with you when calling a poison control center or doctor, or going for treatment.	

100405

**STORAGE AND DISPOSAL**

Do not contaminate water, food, or feed by storage and disposal.

**Pesticide Storage:** Store away from heat.

**Container Disposal:** Do not reuse empty container.

Discard in trash.

EPA REG. NO. 10492-5-8155 EPA EST. NO. 8155-OH-1



This container is coded for Recycling



College of Knowledge

IEHA Certified Training & Education Available on CD ROM



815

**ULTRA SPRAY DISINFECTANT**

**KILLS ALL LISTED PATHOGENS IN ONE MINUTE\***  
**BACTERICIDAL - TUBERCULOCIDAL - VIRUCIDAL\*\***  
**- FUNGICIDAL\*\*\* - HERPECIDAL**

**ACTIVE INGREDIENTS**

N-Alkyl (60% C <sub>12</sub> , 30% C <sub>14</sub> , 5% C <sub>16</sub> , 5% C <sub>18</sub> )	dimethyl benzyl ammonium chloride .....	0.12%
N-Alkyl (68% C <sub>12</sub> , 32% C <sub>14</sub> )	dimethyl ethyl benzyl ammonium chloride ..	0.12%
Isopropyl alcohol .....		63.25%
OTHER INGREDIENTS .....		36.51%
<b>TOTAL</b>		<b>100.00%</b>



**Keep Out of Reach of Children WARNING**

See back panel for additional precautionary statements.

Net Contents: One Quart (32 Fl. Oz.) 946 ml

Manufactured For:  
**CANBERRA CORPORATION**  
 3610 Holland-Sylvania Rd.  
 Toledo, Ohio 43615